

Carrying Cases

Carrying Cases



29-20000
MC400x30
Carrying Case



29-20200
MC1000x9
Carrying Case



29-20300
MC1000x33
Carrying Case



29-20810
Customizable
Foam to fit
MiniCans and
Samplers



29-20400
16 MC1400L
Carrying Case



29-20500
Bottle-Vac™ / MC600
Carrying Case w/
Fitting Compartments



4 Position Bottle-Vac™
Shipping Box(29-21204)

Carrying Cases and Shipping Supplies

Protective dust caps are available for Micro Valves and Miniature Quick Connects to help prevent the introduction of particulates during transportation. Leak-proof caps are available as well to insure very long term storage under vacuum, or when subjecting o-ring sealing valves to cold conditions (below minus 10 deg. C).

Several carrying cases are available for MiniCans, providing options based on the number and size of canisters needing to be shipped. Refer to the chart below for canister and valve compatibility. Custom kits can be made out of the cases with pre-cut foam inserts to accommodate canisters, samplers, and test gauges.

	29-20000	29-20200	29-20300	29-20400	29-20500
MC400L	30		33		15
MC400QT2	30		33		15
MC400V					15
MC600QT2			33		15
MC1000L		9	33		
MC1000V		9	33		
MC1400L				16	
MC1400QT2				16	
MC1400V				16	
BV470QT2					15

Case PN 29-20500 also has compartments for 2 CS1200P Samplers, a vacuum gauge, and 14 40-60ml Vials to hold CS1200 restrictors, restricted samplers, filters, wall samplers and other accessories. The 40-60ml vials can also be used as vacuum samplers when fitted with PN 30-22610 Bottle-Vac™ fittings (see Pg. 20).

Unit	Part No.	Description
100/PK	16-17890	Plastic Miniature Quick Connect Shipping Caps
100/PK	16-17840	Plastic Micro Valve Shipping Caps
EA	30-22070	Micro Valve Leak Tight Cap and Chain
EA	29-20000	Carrying Case for 30 MC400 MiniCans -27.5" x 18" x 8.5"
EA	29-20200	Carrying Case for 9 MC400 Dual End or MC1000 MiniCans
EA	29-20300	Carrying Case for 33 MC1000 MiniCans - 25" x 20" x 18"
EA	29-20400	Carrying Case for 16 MC1400 MiniCans
EA	29-20500	Carrying Case for Bottle-Vacs and Sampling Devices
EA	29-20810	Customizable Foam Case. 16.9"L x 10"W x 6.5"D
EA	29-21104	4 Position, 8oz Bottle-Vac™ Shipping Box with Foam
EA	29-21204	4 Position, 16oz Bottle-Vac™ Shipping Box with Foam