


# Silonite® Coated LOOPS & CUT TUBING



*Silonite® Coated Loops*

Silonite® Coated Loops are pre-bent before coating. This creates the most inert internal surface possible. Loops are batch tested for inertness using dry, polar standards at room temperature.

## Applications

- ◀ Sulfur Compound analysis
- ◀ Formaldehyde analysis
- ◀ Polar Compound analysis
- ◀ Faster cleanup

Part No.	Description
07-60050	50 ul Deactivated Silonite Loop
07-60125	125 ul Deactivated Silonite Loop
07-60250	250ul Deactivated Silonite Loop
07-60500	0.5cc Deactivated Silonite Loop
07-61000	1.0cc Deactivated Silonite Loop
07-62000	2.0cc Deactivated Silonite Loop
07-65000	5.0cc Deactivated Silonite Loop

*Includes nuts and ferrules*



*Silonite® Coated Cut Tubing*

Small lengths of tubing can be coated on both the inside and outside for use in sampling systems and analyzers.

## Features

- ◀ Superior inertness, inside and outside
- ◀ Pre-bent to customer specifications

Part No.	Description	Max. Length
15-88020	1/32" OD x 0.022"	ID 10'
15-88120	1/16" OD x 0.02"	ID 8'
15-88140	1/16" OD x 0.04"	ID 20'
15-88280	1/8" OD x 0.083"	ID 50'