

Silonite[®] Coated INJECTION LINERS




Silonite[®] Coated Injection Liners

Capillary Injection Liners

Silonite[®] coated GC Injection Liners provide a "no-bleed" solution for the analysis of reactive compounds. Available to fit most common GC models.

Available in Packages of 5 and in Bulk.
Contact Entech for Ordering Information:
805-527-5939


Injection Liners for Agilent Instruments



Part No.	Description
Agilent 5890/6890	
37-45010	Split/Splitless FocusLiner™, 78.5 x 6.3mm, 4mm ID
37-45020	Split/Splitless Tapered FocusLiner™, 78.5 x 6.3mm, 4mm ID
37-45030	Split/Splitless Fast FocusLiner™, 78.5 x 6.3mm, 2.3mm ID
37-45040	Split/Splitless Tapered Fast FocusLiner™, 78.5 x 6.3mm, 2.3mm ID
37-45050*	Split Straight-through Liner, 78.5 x 6.3mm, 4mm ID
37-45060*	Split w/Quartz Wool, 78.5 x 6.3mm, 4mm ID
37-45070*	Split/Splitless w/Single Taper, 78.5 x 6.3mm, 4mm ID
37-45080*	Split/Splitless w/Single Taper, Quartz Wool, 78.5 x 6.3mm, 4mm ID
37-45090*	Direct, Straight-through Liner, 78.5 x 6.3mm, 1.2mm ID
37-45100	Splitless (Quartz) Straight-through Liner, 78.5 x 6.3mm, 2mm ID
37-45110*	Splitless w/Recessed Gooseneck, 78.5 x 6.3mm, 2mm ID


* Also fits Varian 1177 Injector

Injection Liners for Shimadzu Instruments




Part No.	Description
17A (SPL-17 Injectors)	
37-43150	Split/Splitless, FocusLiner™, 95.5 x 5.0mm, 3.4mm ID
37-43010	Split/Splitless, Tapered FocusLiner™, 95.5 x 5.0mm, 3.4mm ID
37-43110	Split/Splitless w/Recessed Gooseneck & Quartz Wool, 95.5 x 5.0mm, 3.4mm ID
37-43130	Split, Straight-through Liner, 95.5 x 5.0mm, 3.4mm ID
37-43140	Splitless, Straight-through Liner, 95.5 x 5.0mm, 2.6mm ID
37-43160	Direct, 95.5 x 5.0mm, 2.3mm ID
37-43210	Split/Splitless w/Middle Gooseneck, 95.5 x 5.0mm, 3.4mm ID

Injection Liners for Perkin-Elmer Instruments




Part No.	Description
Perkin-Elmer AutoSystem™	
37-46210	Split/Splitless FocusLiner™, 92 x 6.2mm, 4mm ID
37-46030	Split/Splitless, Tapered FocusLiner™, 92 x 6.2mm, 4mm ID
37-46110	Split, Straight-through Liner, 92 x 5.2mm, 4mm ID

Injection Liners for Thermo Finnigan Instruments



Part No.	Description
Mega Series 4000/5000/8000	
37-44010	Split/Splitless, Recessed Gooseneck w/Wool, 79.5 x 5mm, 3mm ID
37-44020	Splitless, Straight-through Liner, 79.5 x 5mm, 2mm ID
Model 8000/TRACE™	
37-44030	Split/Splitless, FocusLiner™, 105 x 8mm, 5mm ID
37-44040	Split/Splitless, Tapered FocusLiner™, 105 x 8mm, 5mm ID
37-44050	Splitless w/Single Taper, 105 x 8mm, 3mm ID
37-44060	Split w/Single Taper, 105 x 8mm, 5mm ID
37-44070	Split, Straight-through Liner, 105 x 8mm, 5mm ID

Injection Liners for Varian Instruments



Part No.	Description
Injector Models 1075/1077 - GC Model 3400 & 3600	
37-42130	Split/Splitless, FocusLiner™, 72 x 6.3mm, 4mm ID
37-42210	Split/Splitless, Tapered FocusLiner™, 72 x 6.3mm, 4mm ID
37-42310	Split/Splitless, FocusLiner™ w/Top-end Restriction, 72 x 6.3mm, 4mm ID
37-42110	Split/Splitless, Fast FocusLiner™, 72 x 6.3mm, 2.3mm ID
37-42150	Split w/Quartz Wool, 72 x 6.3mm, 4mm ID
Injector Models 1078/1079 - Temperature Programmable Injector	
37-42160	Split/Splitless, FocusLiner™, 54 x 5mm, 3.4mm ID
37-42180	Split/Splitless, Tapered FocusLiner™, 54 x 5mm, 3.4mm ID
37-42200	Split/Splitless w/Single Taper, 54 x 5mm, 3.4mm ID
37-42220	Splitless w/Single Taper, 54 x 5mm, 2mm ID
37-42240	SPI Liner (restriction 0.25mm), 54 x 4.6mm, 0.5mm ID
37-42260	SPI Liner (restriction 0.25mm, for 0.53mm ID on-column), 54 x 4.6mm, 0.5mm ID

Injector Model 1177 & GC Models 3380/3800/3900 - See Agilent Injection Liners

Available in packages of 5 and in bulk.
 Contact Entech for ordering information:
 805-527-5939